


SRVTM DATA SHEET

- Single Use
- Connects using standard Luer Connectors
- Makes return to Treatment fast, simple, & clean
- Proprietary Flow Path Seals when Valve is separated under predetermined tension



SAFETY RELEASE VALVE™ - Performance Specifications	
Luer Connectivity	Yes
Single Use Only	Yes
Materials	Silicone, Polycarbonate
Priming Volume	0.182 mL
Residual Volume following Activation	0.058 mL
Non-Metallic	Yes
Flow Path	In-line proprietary design
Latex/DEHP Free	Yes
Disconnect Tension	Max 3.25 lbs

SAFETY RELEASE VALVE™ - Order Information	
1812-B300 Rev D	
Individual pkgd unit	
100 units	

Rx Only. Copyright 2022 Linear Health Sciences, LLC Oklahoma City, OK All Rights Reserved Ref: SRV-M-2 Rev 03

